

# Overview

Using DXLite software you can send static or dynamic digital radiographs over the internet for digitizing services. All you need is a

1. Windows XP or Vista system,
2. Cable/DSL/Fiber optics,
3. Digital Camera/Scanner to convert static X-rays films to digital format (JPG, JPEG, BMP or TIFF) or
4. Motion Video recorded in internet format (AVI or MPEG).

DXLite is very easy to use and intuitive with

- Wizard based submissions
- Drag 'n' drop images
- Suspend and resume submissions
- Patient cloning
- Scheduled X-ray image or movie uploads
- Scheduled full color report downloads and Status notifications

While it is advanced enough to

- Submit multiple patients in one single session.
- Pool large size movies and schedule upload activity
- Perform automatic updates.
- Customize reports for better presentation and printing performance
- Backup and restore patient database.

No more paper work! No FedEx shipping delays! After downloading and registering your personal copy of the DXLite you can send x-rays over the SSL encrypted, HIPAA compliant website anytime (day or night) and get high quality reports for better treatment.